




Universal Nasal Home Collection Kit Instructions

For use under Emergency Use Authorization



Each kit bag contains:

One (1) transport vial containing media, one (1) swab,
one (1) biosafety bag, and one (1) absorbent pad.

*Kit materials may differ in appearance based on what supplies are
available at a given time.*

1



Watch the Instructional Video

Visit www.unh.edu/covidtest and watch
the instructional video.

2



Label

Fill out one barcode label with the date
and time of sample collection. You
should receive a sheet of personalized
labels when you pick up your kits.

3



Prepare

Affix one label on the vial. Cover the barcode
that is already on the vial.

**You no longer need to place a label on the
biosafety bag.**


4



Prior to Test

Wash hands with soap and
water for 20 seconds and dry,
or use hand sanitizer.


5



Prep for Collection

Remove swab from packaging.
Open from the side opposite the
swab tip.


6



Sampling

Rotate the swab in one
nostril for 15 seconds
(as shown in the video).


7



Sampling

Rotate the swab in
the other nostril for
15 seconds.

7



Swab Securing

Place the swab (tip down) in the vial. Snap the
swab handle by bending it firmly against the lip of
the vial and dispose of the remaining piece. Screw
the lid on the vial to seal securely. Make sure the
top of the handle sits below the lip of the vial so
the handle is not bent when the lid is sealed.
NEVER fold the swab handle over to fit it in the vial.


8



Sample Securing

Place the sealed tube in the
biosafety bag along with the
absorbent pad and seal the bag. (In
some kits, the absorbent pad may
already be inside the biosafety bag.)

9



Sample Drop-Off

Drop the sealed biosafety bag in the secured
drop box at any collection point **within 2
hours of sample collection**. Please use the
drop boxes located on the campus where you
received your kits (e.g., if you received your kits
in Durham, please use the Durham drop boxes).



Test kit drop box locations:

Durham Campus

- > A-Lot tunnel
- > Dimond Library: Under the front covered walkway
- > Hamel Recreation Center: Lobby
- > Health & Wellness: Inside main entrance
- > Holloway Commons (HoCo): By turnstiles inside main door
- > Kingsbury: Albert's Café Lobby
- > MUB: 213 COVID Office
- > Philbrook Dining Hall: By turnstiles inside main door
- > Stillings Dining Hall: Main entrance (locks Friday at noon)
- > West Edge Hut

UNH Manchester Campus

- > Outside Front Entrance of 88 Commercial Street

UNH School of Law Campus

- > Outside Main Lobby of 2 White Street

Additional drop box information:

- > Drop boxes will open by 9 am on Monday and be closed on Friday's at 4 pm. (**UNH Manchester drop box remains open at all times.**)
- > The last run of each day will be 4 pm. Anything left in the box (Monday – Thursday) after 4 pm will be processed the following day.
- > Any samples left on top of or around a locked box (4 pm Friday – 9 am Monday) will be taken to the lab and immediately rejected.
- > The lab will be closed on university holidays and all boxes will reopen by 9 am the next morning.

Visit www.unh.edu/covidtest for more details.

Frequently asked questions

What happens if I test positive for COVID-19?

- > UNH Health & Wellness will contact you with further instructions. They will also contact you daily to monitor your symptoms, either by voice call, text or email.
- > Health & Wellness will help you connect with Academic Affairs to answer questions about classes

What happens if I forget to drop off my self-test sample on time?

- > In order to maintain a valid Wildcat Pass, you must test in accordance with your testing schedule. Failure to do so could result in conduct charges through the Office of Community Standard

How can I report my vaccination status?

- > You can upload a vaccine card at covid-19.unh.edu/UploadVaccine.html

What happens if the test tube leaks when I flip it over to check the seal?

- > You need to complete a new testing sample.

What if my sample spills, is rejected, or I've run out of test?

- > Test kits will be provided periodically throughout the semester, but you will not be given any extra. If your sample has been compromised or a kit is no longer usable, please contact Health & Wellness at (603)-862-9355.

What is the absorbent pad for?

- > The absorbent pad is included to soak up spilled liquid in the event a sample leaks. This helps protect lab workers and prevents the contamination of other samples.

What happens if I am experiencing COVID-19 related symptoms?

- > You should self-isolate immediately and call Health & Wellness at (603) 862-9355 to be evaluated.

How will I eat if I am required to stay in an isolation hall?

- > UNH Dining will prepare boxed meals and deliver them to you. They will take into consideration dietary needs.

Contact information:

Visit UNH Health & Wellness website for more information at www.unh.edu/health or contact covid@unh.edu.

For issues with your testing portal, contact unhcnvlab@unh.edu

For all COVID related academic concerns, contact the Associate Dean of your college.



COVID Incident Reporting Form
tinyurl.com/UNHCovid



Upload your vaccine card
covid-19.unh.edu/UploadVaccine.html